

 Science **made** smarter

Welcome to Interacoustics
**Where
Innovation
Meets
Audiology
Excellence**



Interacoustics

Audiometry

Tympanometry



ABR

OAE

Hearing Aid Fitting

Balance

Contents



Welcome to Interacoustics

Your Leading Provider of Screening and Diagnostic Audiology Equipment

At Interacoustics, we take pride in being at the forefront of audiology innovation, dedicated to providing top-notch screening and diagnostic solutions. With a commitment to excellence, we offer a comprehensive range of cutting-edge equipment designed to meet the evolving needs of healthcare professionals.

Our Commitment to Quality:

- **Precision and Accuracy:** Our audiology equipment are engineered for precision, ensuring reliable results for accurate diagnoses.
- **User-Friendly Design:** Experience the ease of use with our equipment. We understand that simplicity is key, and our tools are crafted to seamlessly integrate into your workflow, making your tasks more efficient and effective.
- **Innovation:** Stay ahead in the field with our state-of-the-art technology, designed to enhance your diagnostic capabilities and patient care.

Explore Our Offerings:

- Discover a wide array of audiology solutions tailored to your specific requirements:
- **Screening Equipment:** Effortlessly integrate efficient screening processes into your workflow.

Experience the Interacoustics Advantage:

- **Reliability:** Trust in the reliability of our equipment, backed by years of expertise in audiology solutions.
- **Customer Support:** Benefit from dedicated customer support, ensuring a seamless experience throughout your partnership with Interacoustics. With a team of 3 audiologists and 7 technical support staff we will have all bases covered with getting you up and running.

For detailed information on our products, services, and how to partner with us, visit our website at www.interacoustics.com.

Contact Us:

We look forward to assisting you on your journey. Feel free to reach out to our team at info@interacoustics.com.au or (02) 8899 1200.

Additional services



Consumables:

Sanibel Supply specialises in supplies for all audiology equipment, clinical consumables, infection control and more. Sanibel's goal is to make it simple and affordable for hearing professionals to buy supplies. Sanibel Supply produces its own line of quality products for most of the major brands in the industry, in addition to stocking consumables for all clinical needs.

- ADI tips
- speculum
- Probe tip
- insert tips
- infection control
- handheld otoscope and tips
- single use flyer/information

[Sanibel online catalogue](#)

Service:

Ensure the Accuracy of Your Clinical Equipment!

At Interacoustics we understand the critical importance of precision in clinical settings. To maintain the highest standards, we recommend annual calibration for all clinical equipment, including Tympanometers and audiometers from Diatec.

Why Annual Calibration?

- **Accuracy Assurance:** Annual calibration ensures that your equipment maintains optimal accuracy, providing reliable results for accurate diagnoses.
- **Compliance:** Stay compliant with industry standards and regulations by scheduling regular calibrations for your clinical instruments.

Booking and Fees Information:

For your convenience, all information regarding bookings and fees can be found on the Diatec website: www.diatec.com.au.

How to Book:

1. Visit the website: www.diatec.com.au
2. Navigate to the "Calibration Services" section.
3. Follow the easy steps to schedule your annual calibration.

Additional Services:

In addition to calibration services, Diatec offers:

- **Noise Assessment Services:** Ensure your testing environment is suitable for accurate hearing assessment.
- **General Repairs of Equipment:** Trust our experts to handle the general repairs of your clinical equipment promptly and efficiently.

For inquiries and more information, feel free to contact Diatec on info@diatec.com.au or (02) 8899 1280

Products Audiometers

AS608 b and AS608 e

Quick hearing screening **made real**



AS608

A portable audiometer made for schools, industry and medical practices; wherever quick hearing screening evaluations are performed.

AS608 b

Specifics:

- Weighs less than 1.6kg
- Manual audiometry
- Battery or power supply operated
- Dedicated carrying bag

Basic AS608

- 11 frequencies from 125 Hz to 8 kHz
- Customize frequency test selections
- Output levels from -10 dB to 100 dB HL
- Pure tone, warble tone and pulsed tone
- Silent stimulus switch
- Patient response switch
- Power supply

AS608 e

Specifics - The same inclusions as the basic unit PLUS:

- USB Cable
- PC transfer
- Automatic HW testing

Additional Features - The same inclusions as the basic unit PLUS:

- Talk Forward
- Hughson Westlake automatic pure tone test
- PC integration with Diagnostic Suite

[Packing list](#)

[Read more](#)

Products Audiometers

Luna

Audiometry **made** plug and play

Perform fast and reliable screening audiometry with the Luna screening audiometer.

Plug in your patient response to your Luna USB headset and connect the headset to your Windows® tablet or laptop and you are ready to start testing your patients



Portable

The headset audiometer only weighs 500 grams including the patient response.



Fast


Test at scale with automated tests and plug-and-play functionality.



Reliable

Randomized frequency order for reliable results and no bulky equipment to manage.

[Read more](#)



Products Audiometers

AD226

Air and bone audiometry **made** easy

AD226 is made to perform routine air and bone conduction threshold audiometry. For advanced audiometry, you can add several hearing tests to the test battery with the extended license.

- Storage space for 500 patients
- Portable: option to run on AA batteries
- Lightweight with optional carrying bag

Option:

- Option to use via PC



[Read more](#)

Products Audiometers

PA 5

Paediatric testing **made** comfortable

The PA5 is a free field audiometer designed for trained healthcare professionals involved in testing young children. It features one-handed operation, silent buttons and a variety of light and sound stimuli.

- Conditioned Orientation Response (COR)
- Visual Reinforcement Audiometry (VRA)
- Auropalpebral Reflex (APR)

Features:

- Flashing Light Stimulation
- TDH39 Audiometric Headphone (Optional)
- Robust, Lightweight
- Standard AA Batteries
- Automatic Battery Switch
- Automatic Battery Indication



[Read more](#)

Products Tympanometers

MT10

Tympanometry **made** portable

MT10b

The MT10b is our entry level tympanometer which performs 226Hz tympanometry only. It comes supplied with:

- Batteries
- An initial box of disposable probe tips, and
- Cleaning accessories

MT10 e

The MT10e is our entry level tympanometer which performs 226Hz tympanometry and ipsilateral reflexes. It comes supplied with:

- Batteries
- An initial box of disposable probe tips, and
- Cleaning accessories.

Optional Thermal Printers or Infrared Dongles for PC upload are available with the MT10 e model.



[Packing list](#)

[Read more](#)

Products Tympanometers

Titan Screener

Titan delivers tympanometry through advanced technologies. These include: "Endless airflow" pressure control technology, high resolution recordings and an intelligent pump system with an adaptive speed control around the tympanic peak.

Titan screener features:

- 226, 678, 800 & 1000Hz tympanometry
- Ipsilateral reflexes
- PC transfer of results
- Titan Suite Software



[Packing list](#)

[Read more](#)

Products Tympanometers - Comparison Table



Titan Screener

MT10E

	Titan Screener	MT10E
226HZ	✓	✓
IPSI Reflexes	✓	✓
HF Tymps	✓	
Easy PC upload	✓	
Ability to add on OAE	✓	
Price point	\$\$	\$
Carry Case	✓	Optional

Products Video Oscope

Video Oscope

Video otoscopy **made** effective

The Viot™ from Interacoustics is the natural choice for clinics that need an easy-to-use video otoscope for preliminary diagnostic assessments. Viot™ objectively views the external auditory canal before, during and after any audiological procedure.

A video otoscope with high usability

Weighing in at 135 g, Viot™ is truly a compact and hand-held video otoscope. The built-in LED provides a powerful and bright light source, allowing true color capturing. Finally, you do not need an additional power supply, as the video otoscope is powered via its USB cable, which is a USB **2.0 connection**.



[Read more](#)

Training section

Interacoustics prides itself on the educational resources available for our customers. Our product and Academy resources are world leading and we encourage you to browse the available training videos and documents.

Blogs: <https://interacoustics.com.au/updates/blog>

Packing Lists: <https://interacoustics.com.au/resources>

YouTube Channel: <https://www.youtube.com/user/InteracousticsTV>

Academy training: <https://www.interacoustics.com/academy>

Interacoustics Australia also offers tailored training programs for companies seeking personalised service. Please contact us at info@interacoustics.com.au to discuss a suitable program relevant to your needs.

We regularly provide online education courses and participate in conferences. To keep up to date with all our latest news and events, please sign up to our newsletter on our website home page: www.interacoustics.com.au.

T&C

Warranty

Returns etc.

Warranty & Guarantee

New and unused goods delivered directly by the supplier or its authorised representatives carry a guarantee under the following conditions.

- When properly installed and used the goods are guaranteed for two (2) years from the date of delivery. This means that the supplier will correct, without further charge, defects in the goods arising from faulty manufacture or deficient materials or components, providing the goods are returned, postage paid to the supplier or the supplier's authorised representative before the expiry of the guarantee period together with a full account of the defects.
- The cost of packaging and shipping the repaired goods back to the purchaser are covered by the supplier or the supplier's authorised representative.
- Parts and accessories are warranted for ninety (90) days.
- The supplier shall not be held liable for any loss or injury incurred as a result of repairs or modifications made to the goods by anyone other than the supplier or the supplier's authorised representative without the supplier's written permission. Such repairs or modifications will also terminate any further guaranteed obligations.
- Compensation cannot be paid for wages, loss of profit, waste of time, packing and transport resulting from defects in the goods.

Payments

Interacoustics will issue an invoice for payment. Payments are to be made thirty (30) days from the date of invoice unless otherwise agreed in writing with a sales representative.

If payment is not made within the payment terms interest may be charged on the unpaid amount.

Delivery and Returns

Rates

Our goods are shipped via air freight from manufacturing plants in Europe. Most items we provide are custom produced. This allows the flexibility to design your own unit to suit your clinical needs. Our delivery time frames will be discussed with you at the time of purchase. All reasonable efforts will be made to give you a clear indication of the wait time. Recent supply chain and production issues as well as local freightage and courier delays have pushed out delivery timeframes. Some items may need to be occasionally sea freighted.

Unless an item is in stock in our local storage facilities, units take a minimum of 4 - 6 weeks to be delivered from your point of order.

Freight costs are not included in the price of your instrument and will be quoted at the time of the purchase. Charges are based on volume and delivery location. We will deliver your goods to the nominated address and place. Please provide clear delivery instructions at the time of purchase.

Risk of damage to or loss of the goods shall pass to you at the time of delivery.

Returns

All sales are final and we do not accept returns or offer refunds for a change of mind. On an individual basis, items may be returned at the sole discretion of the General Manager. A restocking fee of 15-20% applies based on the item purchased.

The rights given under the product warranty are in addition to other rights and remedies under a law in relation to the goods to which the warranty relates to. Our goods come with guarantees that cannot be excluded under the Australian Consumer Laws. You may be entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You may also be entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Received a wrong item?

We do all we can to deliver your products in the correct quantity and condition. If you find that we have made a mistake, we kindly ask you to send an e-mail to info@interacoustics.com.au. A customer service representative will contact you to correct the problem.

Contacting Us

If at any time you would like to contact us to discuss our Terms and Conditions, you can do so by sending an e-mail to us at info@interacoustics.com.au

Science made smarter

Interacoustics is more than state-of-the-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

Interacoustics.com

Interacoustics Pty Ltd

Level 5, 11 Khartoum Road

Macquarie Park NSW 2113

Ph: +612 8899 1200

info@interacoustics.com.au

